



Mouse anti-Human MSL2 monoclonal antibody, clone 4B5 (CABT-B9476)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	MSL2 (NP_060603, 478 a.a. ~ 575 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	4B5
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords

MSL2; male-specific lethal 2 homolog (Drosophila); MSL-2; MSL2L1; RNF184; E3 ubiquitin-protein ligase MSL2; MSL2-like 1; ring finger protein 184; male-specific lethal 2-like 1; male-specific lethal-2 homolog 1;

GENE INFORMATION

Entrez Gene ID[55167](#)**UniProt ID**[Q9HCL7](#)