



Mouse anti-Human NACA2 monoclonal antibody, clone 5I7 (CABT-B9482)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	NACA2 (NP_954984, 116 a.a. ~ 215 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	5I7
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords

NACA2; nascent polypeptide-associated complex alpha subunit 2; ANAC; NACAL; nascent polypeptide-associated complex subunit alpha-2; alpha-NAC protein; hom s 2.01 protein; IgE autoantigen alpha-nascent polypeptide-associated complex; nascent polypeptide-associated complex, alpha polypeptide, 2;

GENE INFORMATION

Entrez Gene ID

[342538](#)

UniProt ID

[Q9H009](#)
