



Mouse anti-Human NCS1 monoclonal antibody, clone 4F8 (CABT-B9484)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	NCS1 (AAH04856, 1 a.a. ~ 190 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	4F8
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords

NCS1; neuronal calcium sensor 1; Freq; Mfreq; NCS-1; AI836659; 9430075O15Rik;
A730032G13Rik; frequenin homolog; neuronal calcium sensor-1;

GENE INFORMATION

Entrez Gene ID[14299](#)**UniProt ID**[Q8BNY6](#)