



# Mouse anti-Human NDFIP2 monoclonal antibody, clone 3D8 (CABT-B9485)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	NDFIP2 (AAH21988, 1 a.a. ~ 242 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
<b>Isotype</b>	IgG2a, κ
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	3D8
<b>Purification</b>	Protein A
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, WB
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	PBS, pH 7.4
<b>Preservative</b>	no preservative
<b>Storage</b>	-20°C, Avoid Freeze/Thaw Cycles

## BACKGROUND

**Keywords** NDFIP2; Nedd4 family interacting protein 2; N4WBP5A; NEDD4 family-interacting protein 2;

NF-kappa-B-activating protein 413; NEDD4 WW domain-binding protein 5A; putative NF-kappa-B-activating protein 413; MAPK-activating protein PM04 PM05 PM06 PM07; putative MAPK-activating protein PM04/PM05/PM06/PM07;

---

## GENE INFORMATION

Entrez Gene ID	<a href="#">54602</a>
----------------	-----------------------

---

UniProt ID	<a href="#">B4DGY6</a>
------------	------------------------

---