



# Mouse anti-Human NME5 monoclonal antibody, clone 5G6 (CABT-B9486)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	NME5 (NP_003542, 121 a.a. ~ 212 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
<b>Isotype</b>	IgG2a, κ
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	5G6
<b>Purification</b>	Affinity chromatography
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, IM, WB
<b>Format</b>	Liquid
<b>Concentration</b>	0.2-1.0 mg/ml
<b>Size</b>	100 µg
<b>Buffer</b>	PBS, pH 7.4
<b>Preservative</b>	no preservative
<b>Storage</b>	-20°C, Avoid Freeze/Thaw Cycles

## BACKGROUND

**Keywords**

NME5; NME/NM23 family member 5; NM23H5; RSPH23; NM23-H5; nucleoside diphosphate kinase homolog 5; NDK-H 5; IPIA-beta; NDP kinase homolog 5; radial spoke 23 homolog; testis-specific nm23 homolog; inhibitor of p53-induced apoptosis-beta; non-metastatic cells 5, protein expressed in (nucleoside-diphosphate kinase);

---

## GENE INFORMATION

**Entrez Gene ID**[8382](#)**UniProt ID**[P56597](#)