



Mouse anti-Human NOL4 monoclonal antibody, clone 3B20 (CABT-B9487)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	NOL4 (NP_003778, 425 a.a. ~ 524 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	3B20
Purification	Protein A
Conjugate	Unconjugated
Applications	ICC, IF, IHC-P
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords NOL4; nucleolar protein 4; NOLP; CT125; HRIHFB2255; cancer/testis antigen 125; nucleolar

localized protein;

GENE INFORMATION

Entrez Gene ID	8715
----------------	----------------------

UniProt ID	B4DLW2
------------	------------------------
