



Mouse anti-Human NUBP1 monoclonal antibody, clone 3C22 (CABT-B9488)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	NUBP1 (NP_002475.1, 222 a.a. ~ 320 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	3C22
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction	NUBP1 is a member of the NUBP/MRP subfamily of ATP-binding proteins (Nakashima et al., 1999 [PubMed 10486206]).[supplied by OMIM, Mar 2008]
Keywords	NUBP1; nucleotide binding protein 1; NBP; NBP1; NBP35; cytosolic Fe-S cluster assembly factor NUBP1; NBP 1; nucleotide-binding protein 1; nucleotide binding protein (e.coli MinD like); nucleotide binding protein 1 (E.coli MinD like); nucleotide binding protein 1 (MinD homolog, E. coli);

GENE INFORMATION

Entrez Gene ID	4682
UniProt ID	P53384