



Mouse anti-Human NXPH3 monoclonal antibody, clone 5D9 (CABT-B9489)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	NXPH3 (AAH22541, 1 a.a. ~ 252 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	5D9
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords

NXPH3; neurexophilin 3; NPH3; neurexophilin-3;

GENE INFORMATION

Entrez Gene ID[11248](#)

UniProt ID[O95157](#)
