



Mouse anti-Human ODAM monoclonal antibody, clone 3G9 (CABT-B9490)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	ODAM (AAH17796.1, 1 a.a. ~ 153 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	3G9
Purification	Protein A
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords ODAM; odontogenic, ameloblast associated; APIN; odontogenic ameloblast-associated

protein;

GENE INFORMATION

Entrez Gene ID	54959
----------------	-----------------------

UniProt ID	A1E959
------------	------------------------
