



Mouse anti-Human PAPOLB monoclonal antibody, clone 4E8 (CABT-B9497)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	PAPOLB (AAH36653, 1 a.a. ~ 637 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	4E8
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, ICC, IF, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords PAPOLB; poly(A) polymerase beta (testis specific); PAPT; TPAP; poly(A) polymerase beta; PAP-beta; testis-specific poly(A) polymerase; polynucleotide adenylyltransferase beta;

GENE INFORMATION

Entrez Gene ID [56903](#)

UniProt ID [A4D1Z6](#)
