



Mouse anti-Human PDIK1L monoclonal antibody, clone 5C8 (CABT-B9498)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	PDIK1L (NP_690048, 1 a.a. ~ 99 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	5C8
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords

PDIK1L; PDLIM1 interacting kinase 1 like; CLIK1L; STK35L2; serine/threonine-protein kinase
PDIK1L; casein kinase; PDLIM1-interacting kinase 1-like;

GENE INFORMATION

Entrez Gene ID[149420](#)**UniProt ID**[Q8N165](#)