



Mouse anti-Human PIPOX monoclonal antibody, clone 4E2 (CABT-B9502)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	PIPOX (AAH27622.1, 292 a.a. ~ 384 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	4E2
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, ICC, IF, IM, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords PIPOX; pipecolic acid oxidase; LPIPOX; peroxisomal sarcosine oxidase; PSO; L-pipecolate oxidase; L-pipecolic acid oxidase;

GENE INFORMATION

Entrez Gene ID [51268](#)

UniProt ID [Q9P0Z9](#)
