



Mouse anti-Human RASL12 monoclonal antibody, clone 2D20 (CABT-B9515)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	RASL12 (AAH53734, 1 a.a. ~ 266 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2b, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	2D20
Purification	Protein A
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords RASL12; RAS-like, family 12; RIS; ras-like protein family member 12; ras-like protein Ris; Ras

family member Ris;

GENE INFORMATION

Entrez Gene ID [51285](#)

UniProt ID [A0A024R5Y3](#)
