



# Mouse anti-Human RASSF9 monoclonal antibody, clone 3G9 (CABT-B9516)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	RASSF9 (AAH31589, 1 a.a. ~ 435 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	3G9
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

## BACKGROUND

**Introduction**

The protein encoded by this gene localizes to perinuclear endosomes. This protein associates with peptidylglycine alpha-amidating monooxygenase, and may be involved with the trafficking of this enzyme through secretory or endosomal pathways. [provided by RefSeq, Jul 2008]

**Keywords**

RASSF9; Ras association (RalGDS/AF-6) domain family (N-terminal) member 9; PAMCI; PCIP1; P-CIP1; ras association domain-containing protein 9; PAM COOH-terminal interactor protein 1; peptidylglycine alpha-amidating monooxygenase COOH-terminal interactor protein-1;

## GENE INFORMATION

**Entrez Gene ID**

[9182](#)

**UniProt ID**

[O75901](#)