



Mouse anti-Human RhoBTB3 monoclonal antibody, clone 4B20 (CABT-B9523)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	RHOBTB3 (NP_055714.2, 1 a.a. ~ 99 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	4B20
Purification	Protein A
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction RHOBTB3 is a member of the evolutionarily conserved RHOBTB subfamily of Rho GTPases.

For background information on RHOBTBs, see RHOBTB1 (MIM 607351).[supplied by OMIM, Apr 2004]

Keywords	RHOBTB3; Rho-related BTB domain containing 3; rho-related BTB domain-containing protein 3;
-----------------	--

GENE INFORMATION

Entrez Gene ID	22836
-----------------------	-----------------------

UniProt ID	O94955
-------------------	------------------------
