



Mouse anti-Human RHOXF1 monoclonal antibody, clone 5E23 (CABT-B9524)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	RHOXF1 (NP_644811.1, 1 a.a. ~ 110 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	5E23
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, ICC, IF, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction

This gene is a member of the PEPP subfamily of paired-like homeobox genes. The gene may be regulated by androgens and epigenetic mechanisms. The encoded nuclear protein is likely a transcription factor that may play a role in human reproduction. [provided by RefSeq, Dec 2012]

Keywords

RHOXF1; Rhox homeobox family, member 1; OTEX; PEPP1; rhox homeobox family member 1; PEPP subfamily gene 1; paired-like homeobox protein PEPP-1; ovary-, testis- and epididymis-expressed gene protein;

GENE INFORMATION

Entrez Gene ID

[158800](#)

UniProt ID

[Q8NHV9](#)