



Mouse anti-Human RHOXF2 monoclonal antibody, clone 4D5 (CABT-B9525)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	RHOXF2 (AAH21719.1, 1 a.a. ~ 288 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	4D5
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords

RHOXF2; Rhox homeobox family, member 2; THG1; CT107; PEPP2; PEPP-2; rhox homeobox family member 2; PEPP subfamily gene 2; testis homeobox gene 1; cancer/testis antigen 107; homeobox protein from AL590526; paired-like homeobox protein PEPP-2;

GENE INFORMATION

Entrez Gene ID[84528](#)**UniProt ID**[Q9BQY4](#)