



Mouse anti-Human RNF121 monoclonal antibody, clone 3H8 (CABT-B9527)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	RNF121 (NP_060790.2, 193 a.a. ~ 301 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	3H8
Purification	Protein A
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction The protein encoded by this gene contains a RING finger, a motif present in a variety of

functionally distinct proteins and known to be involved in protein-protein and protein-DNA interactions. Several alternatively spliced transcript variants have been noted for this gene, however, not all are likely to encode viable protein products. [provided by RefSeq, Sep 2008]

Keywords	RNF121; ring finger protein 121; RING finger protein 121;
-----------------	---

GENE INFORMATION

Entrez Gene ID	55298
-----------------------	-----------------------

UniProt ID	C9JQY5
-------------------	------------------------
