



Mouse anti-Human RNF148 monoclonal antibody, clone 3E4 (CABT-B9528)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	RNF148 (NP_932351.1, 139 a.a. ~ 238 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2b, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	3E4
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords RNF148; ring finger protein 148; RING finger protein 148;

GENE INFORMATION

Entrez Gene ID [378925](#)

UniProt ID [Q8N7C7](#)
