



Mouse anti-Human RXFP3 monoclonal antibody, clone 2G20 (CABT-B9535)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	RXFP3 (NP_057652, 1 a.a. ~ 80 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	2G20
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords

RXFP3; relaxin/insulin-like family peptide receptor 3; SALPR; RLN3R1; RXFPR3; GPCR135; relaxin-3 receptor 1; RLN3 receptor 1; relaxin 3 receptor 1; G-protein coupled receptor SALPR; g protein-coupled receptor SALPR; relaxin family peptide receptor 3; G-protein coupled receptor GPCR135; somatostatin and angiotensin-like peptide receptor; somatostatin- and angiotensin-like peptide receptor;

GENE INFORMATION

Entrez Gene ID[51289](#)**UniProt ID**[Q9NSD7](#)
