



Mouse anti-Human SETD3 monoclonal antibody (CABT-B9539)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	2-498 of 594 of SETD3 Protein.
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	IM
Format	Liquid
Concentration	0.5 mg/ml
Size	100 µg
Buffer	PBS
Preservative	0.09% sodium azide
Storage	Maintain refrigerated at 2-8°C for up to 1 month. For long term storage store at -20°C

BACKGROUND

Keywords	SETD3; SET domain containing 3; C14orf154; histone-lysine N-methyltransferase setd3; SET domain-containing protein 3;
-----------------	---

GENE INFORMATION

Entrez Gene ID [84193](#)

UniProt ID [A0A024R6K3](#)
