



Mouse anti-Human SNX20 monoclonal antibody, clone 5F0 (CABT-B9550)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	SNX20 (NP_699168.1, 47 a.a. ~ 128 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	5F0
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, ICC, IF
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction	SNX20 interacts with the cytoplasmic domain of PSGL1 (SELPLG; MIM 600738) and cycles PSGL1 into endosomes.[supplied by OMIM, Feb 2010]
Keywords	SNX20; sorting nexin 20; SLIC1; sorting nexin-20; selectin ligand interactor cytoplasmic 1; selectin ligand interactor cytoplasmic-1; selectin ligand-interactor cytoplasmic 1;

GENE INFORMATION

Entrez Gene ID	124460
UniProt ID	Q7Z614