



Mouse anti-Human SNX24 monoclonal antibody, clone 3B22 (CABT-B9551)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	SNX24 (AAH10886, 1 a.a. ~ 169 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	3B22
Purification	Protein A
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords SNX24; sorting nexin 24; SBBI31; PRO1284; sorting nexin-24; sorting nexin 24;

GENE INFORMATION

Entrez Gene ID [28966](#)

UniProt ID [Q9Y343](#)
