



Mouse anti-Human SPANXN3 monoclonal antibody, clone 2G22 (CABT-B9554)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	SPANXN3 (NP_001009609.1, 1 a.a. ~ 141 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	2G22
Purification	Protein A
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords SPANXN3; SPANX family, member N3; CT11.8; SPANX-N3; sperm protein associated with the

nucleus on the X chromosome N3; nuclear-associated protein SPAN-Xn3; cancer/testis antigen family 11, member 8;

GENE INFORMATION

Entrez Gene ID [139067](#)

UniProt ID [Q5MJ09](#)
