



Mouse anti-Human SPRR4 monoclonal antibody, clone 2H6 (CABT-B9558)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	SPRR4 (NP_775103.1, 1 a.a. ~ 79 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	2H6
Purification	Protein A
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords SPRR4; small proline-rich protein 4; small proline rich protein 4;

GENE INFORMATION

Entrez Gene ID [163778](#)

UniProt ID [Q96PI1](#)
