



Mouse anti-Human SYCE1 monoclonal antibody, clone 7H8 (CABT-B9560)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	SYCE1 (NP_963858.1, 1 a.a. ~ 190 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	7H8
Purification	Protein A
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords SYCE1; synaptonemal complex central element protein 1; CT76; C10orf94; cancer/testis

antigen 76;

GENE INFORMATION

Entrez Gene ID	93426
----------------	-----------------------

UniProt ID	Q8N0S2
------------	------------------------
