



# Rat anti-Mouse TDRD1 monoclonal antibody, clone 840307 (CABT-B9562)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	E. coli-derived recombinant mouse TDRD1 Asn831-Lys960
Isotype	IgG2a
Source/Host	Rat
Species Reactivity	Mouse
Clone	840307
Purification	Protein A/G
Conjugate	Unconjugated
Applications	IHC-Fr
Format	Lyophilized
Concentration	0.2 mg/ml
Size	100 µg
Buffer	PBS with 5% trehalose
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

## BACKGROUND

**Introduction**

This gene is similar to a mouse gene that encodes a tudor domain protein. Alternatively spliced transcript variants have been described but their full length sequences have not been determined. [provided by RefSeq, Jul 2008]

**Keywords**

TDRD1; tudor domain containing 1; CT41.1; tudor domain-containing protein 1; cancer/testis antigen 41.1; tudor domain containing protein 1;

## GENE INFORMATION

**Entrez Gene ID**

[56165](#)

**UniProt ID**

[Q9BXT4](#)