



# Mouse anti-Human TDRD3 monoclonal antibody, clone 92660 (CABT-B9563)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                    |                                                                                        |
|--------------------|----------------------------------------------------------------------------------------|
| Immunogen          | 554-616 of 651 of TDRD3 Protein.                                                       |
| Isotype            | IgG                                                                                    |
| Source/Host        | Mouse                                                                                  |
| Species Reactivity | Human                                                                                  |
| Clone              | 92660                                                                                  |
| Purification       | purified                                                                               |
| Conjugate          | Unconjugated                                                                           |
| Applications       | IM                                                                                     |
| Format             | Liquid                                                                                 |
| Concentration      | 0.5 mg/ml                                                                              |
| Size               | 100 µg                                                                                 |
| Buffer             | PBS                                                                                    |
| Preservative       | 0.09% sodium azide                                                                     |
| Storage            | Maintain refrigerated at 2-8°C for up to 1 month. For long term storage store at -20°C |

## BACKGROUND

**Keywords**

TDRD3; tudor domain containing 3; tudor domain-containing protein 3;

---

## GENE INFORMATION

**Entrez Gene ID**[81550](#)**UniProt ID**[Q9H7E2](#)

---