



Rat anti-Mouse TEX11 monoclonal antibody, clone 347845 (CABT-B9564)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	E. coli-derived recombinant mouse ZIP4H/Tex11 Met768-Leu947
Isotype	IgG2b
Source/Host	Rat
Species Reactivity	Mouse
Clone	347845
Purification	Protein A/G
Conjugate	Unconjugated
Applications	WB
Format	Lyophilized
Concentration	0.5 mg/ml
Size	100 µg
Buffer	PBS with 5% trehalose
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction

This gene is X-linked and is expressed in only male germ cells. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Keywords

TEX11; testis expressed 11; TGC1; TSGA3; testis-expressed sequence 11 protein; testis expressed sequence 11;

GENE INFORMATION

Entrez Gene ID

[56159](#)

UniProt ID

[Q8IYF3](#)