



Mouse anti-Human THG1L monoclonal antibody, clone 5G22 (CABT-B9566)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	THG1L (AAH01523, 1 a.a. ~ 269 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	5G22
Purification	Protein A
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords THG1L; tRNA-histidine guanylyltransferase 1-like (*S. cerevisiae*); ICF45; IHG-1; probable

tRNA(His) guanylyltransferase; induced by high glucose-1; interphase cytoplasmic foci protein 45;

GENE INFORMATION

Entrez Gene ID [54974](#)

UniProt ID [Q9NWX6](#)
