



Mouse anti-Human TLE3 monoclonal antibody, clone DM4684 (CABT-B9567)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to Human TLE3
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	DM4684
Purification	Protein A
Conjugate	Unconjugated
Applications	IHC, WB
Format	Liquid
Size	100 µl
Buffer	PBS, pH 7.2, with 40% glycerol
Preservative	0.02% sodium azide
Storage	Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

BACKGROUND

Introduction	This gene encodes a transcriptional co-repressor protein that belongs to the transducin-like enhancer family of proteins. The members of this family function in the Notch signaling pathway
--------------	--

that regulates determination of cell fate during development. Expression of this gene has been associated with a favorable outcome to chemotherapy with taxanes for ovarian carcinoma. Alternate splicing results in multiple transcript variants. Additional alternatively spliced transcript variants of this gene have been described, but their full-length nature is not known. [provided by RefSeq, Sep 2013]

Keywords

TLE3; transducin-like enhancer of split 3; ESG; ESG3; GRG3; HsT18976; transducin-like enhancer protein 3; enhancer of split groucho 3; enhancer of split groucho-like protein 3; transducin-like enhancer of split 3 (E(sp1) homolog, Drosophila); transducin-like enhancer of split 3, homolog of Drosophila E(sp1);

GENE INFORMATION

Entrez Gene ID [7090](#)

UniProt ID [Q04726](#)
