



Mouse anti-Human TRAPPC5 monoclonal antibody, clone 4E0 (CABT-B9569)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	TRAPPC5 (AAH42161, 1 a.a. ~ 188 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	4E0
Purification	Protein A
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords TRAPPC5; trafficking protein particle complex 5; TRS31; trafficking protein particle complex

subunit 5;

GENE INFORMATION

Entrez Gene ID [126003](#)

UniProt ID [Q8IUR0](#)
