



Mouse anti-Human TRIM46 monoclonal antibody, clone 4H22 (CABT-B9573)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	TRIM46 (NP_079334, 451 a.a. ~ 560 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	4H22
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, ICC, IF
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords TRIM46; tripartite motif containing 46; GENEY; TRIFIC; tripartite motif-containing protein 46; gene Y protein; tripartite, fibronectin type-III and C-terminal SPRY motif protein;

GENE INFORMATION

Entrez Gene ID [80128](#)

UniProt ID [Q7Z4K8](#)
