



Mouse anti-Human TRIM47 monoclonal antibody, clone 4D9 (CABT-B9574)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	TRIM47 (NP_258411, 539 a.a. ~ 638 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	4D9
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, ICC, IF, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords TRIM47; tripartite motif containing 47; GOA; RNF100; tripartite motif-containing protein 47; RING finger protein 100; tripartite motif-containing 47; gene overexpressed in astrocytoma protein;

GENE INFORMATION

Entrez Gene ID [91107](#)

UniProt ID [Q96LD4](#)
