



Mouse anti-Human TWF2 monoclonal antibody, clone 3C6 (CABT-B9577)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	PTK9L (AAH00327, 131 a.a. ~ 230 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	3C6
Purification	Protein A
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction The protein encoded by this gene was identified by its interaction with the catalytic domain of

protein kinase C-zeta. The encoded protein contains an actin-binding site and an ATP-binding site. It is most closely related to twinfilin (PTK9), a conserved actin monomer-binding protein. [provided by RefSeq, Jul 2008]

Keywords

TWF2; twinfilin actin-binding protein 2; A6r; A6RP; PTK9L; MSTP011; twinfilin-2; A6-related protein; twinfilin-1-like protein; twinfilin, actin-binding protein, homolog 2; protein tyrosine kinase 9-like (A6-related protein); PTK9L protein tyrosine kinase 9-like (A6-related protein);

GENE INFORMATION

Entrez Gene ID

[11344](#)

UniProt ID

[Q6IBS0](#)
