



Mouse anti-Human UBE2U monoclonal antibody, clone 4E5 (CABT-B9578)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	UBE2U (NP_542783, 9 a.a. ~ 102 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	4E5
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, ICC, IF
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords

UBE2U; ubiquitin-conjugating enzyme E2U (putative); ubiquitin-conjugating enzyme E2 U;
ubiquitin-protein ligase U; ubiquitin carrier protein U;

GENE INFORMATION

Entrez Gene ID[148581](#)**UniProt ID**[Q5VVX9](#)