



# Mouse anti-Human UNC5CL monoclonal antibody, clone 4H8 (CABT-B9579)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	UNC5CL (AAH35284, 1 a.a. ~ 518 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1, $\kappa$
Source/Host	Mouse
Species Reactivity	Human
Clone	4H8
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, ICC, IF, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 $\mu$ g
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

## BACKGROUND

**Keywords**

UNC5CL; unc-5 homolog C (C. elegans)-like; ZUD; MUXA; UNC5C-like protein; protein unc-5 homolog C-like; ZU5 and death domain containing;

---

## GENE INFORMATION

**Entrez Gene ID**

[222643](#)

---

**UniProt ID**

[H8YHX0](#)

---