



Mouse anti-Human WDR41 monoclonal antibody, clone 2G4 (CABT-B9581)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	WDR41 (AAH40241, 1 a.a. ~ 459 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2b, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	2G4
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords

WDR41; WD repeat domain 41; MSTP048; WD repeat-containing protein 41;

GENE INFORMATION

Entrez Gene ID[55255](#)**UniProt ID**[Q9HAD4](#)
