



Mouse anti-Human ZBTB25 monoclonal antibody, clone 3F9 (CABT-B9584)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	ZBTB25 (NP_008908.2, 331 a.a. ~ 435 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	3F9
Purification	Protein A
Conjugate	Unconjugated
Applications	ELISA, ICC, IF
Format	Liquid
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords ZBTB25; zinc finger and BTB domain containing 25; KUP; ZNF46; C14orf51; zinc finger and

GENE INFORMATION

Entrez Gene ID [7597](#)

UniProt ID [P24278](#)
