



Mouse anti-Human ZFAND3 monoclonal antibody, clone 3E3 (CABT-B9586)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	ZFAND3 (NP_068762.1, 2 a.a. ~ 60 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	3E3
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords

ZFAND3; zinc finger, AN1-type domain 3; TEX27; AN1-type zinc finger protein 3; testis expressed sequence 27;

GENE INFORMATION

Entrez Gene ID[60685](#)**UniProt ID**[A0A024RD12](#)