



Mouse anti-Human ZKSCAN4 monoclonal antibody, clone 2E3 (CABT-B9589)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	ZKSCAN4 (AAH14031, 1 a.a. ~ 545 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Isotype	IgG2a, κ
Source/Host	Mouse
Species Reactivity	Human
Clone	2E3
Purification	Affinity chromatography
Conjugate	Unconjugated
Applications	ELISA, ICC, IF, WB
Format	Liquid
Concentration	0.2-1.0 mg/ml
Size	100 µg
Buffer	PBS, pH 7.4
Preservative	no preservative
Storage	-20°C, Avoid Freeze/Thaw Cycles

BACKGROUND

Keywords

ZKSCAN4; zinc finger with KRAB and SCAN domains 4; ZNF307; ZNF427; ZSCAN36;
P1P373C6; zinc finger protein with KRAB and SCAN domains 4; p373c6.1 (novel C2H2 type
zinc finger protein); zinc finger protein 307; zinc finger protein 427;

GENE INFORMATION

Entrez Gene ID[387032](#)**UniProt ID**[Q969J2](#)