



Rabbit anti-Human LZTS3 Polyclonal antibody (CABT-B9601)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: QEFAMKSVGTRTGGGGSQGSFPGPRGSGSGASRERPGRYPSEDKGLANSLYLNGLRGSD HTDVCGNVVGSSGGSSSSGGSDKAPPQYREPS
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Concentration	0.2 mg/ml
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C.

BACKGROUND

Keywords	LZTS3; leucine zipper, putative tumor suppressor family member 3; PROSAPIP1; leucine zipper putative tumor suppressor 3; proSAP-interacting protein 1; ProSAP/Shank-interacting
-----------------	---

GENE INFORMATION

Entrez Gene ID [9762](#)

UniProt ID [O60299](#)
