



Rabbit anti-Human MORN3 Polyclonal antibody (CABT-B9603)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: SLPDQQTGKCRRVYSGWWKGDKKSGYGIQFFGPKEYYEGDWCGSQRSGWGRMYYSNGDIY EGQWENDKPNGEGLRL
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	IHC, WB
Format	Liquid
Concentration	0.6 mg/ml
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C.

BACKGROUND

Keywords	MORN3; MORN repeat containing 3; MORN repeat-containing protein 3; membrane occupation and recognition nexus repeat containing protein;
-----------------	---

GENE INFORMATION

Entrez Gene ID [283385](#)

UniProt ID [Q6PF18](#)
