



Rabbit anti-Human HYKK Polyclonal antibody (CABT-B9607)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: ALSKPTFSEEQASALVESVFGLK/SKVRPLPSYDDQNFHVYVSKTKDGPTEYVLKISNTK ASKNPDLIEVQNHIMFLKA
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Concentration	0.05 mg/ml
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C.

BACKGROUND

Keywords HYKK; hydroxylysine kinase; AGPHD1; 5-hydroxylysine kinase; 5-hydroxy-L-lysine kinase; aminoglycoside phosphotransferase domain containing 1; aminoglycoside phosphotransferase

domain-containing protein 1;

GENE INFORMATION

Entrez Gene ID [123688](#)

UniProt ID [A2RU49](#)
