



Rabbit anti-Human KIAA1257 Polyclonal antibody (CABT-B9611)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: IQWSVSRPTISLAGASMMEIKELIESESLSSLTNILDRQRSQIKGKDSEGRRKIQRRHK KPLAEEEEADPTLTGPRKQSAFSIQLAVMPLLAGTHCLPCSQ
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	IHC, WB
Format	Liquid
Concentration	0.1 mg/ml
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C.

BACKGROUND

Keywords KIAA1257; uncharacterized protein KIAA1257;

GENE INFORMATION

Entrez Gene ID [57501](#)

UniProt ID [Q9ULG3](#)
