



Rabbit anti-Human FAM21A Polyclonal antibody (CABT-B9615)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: VRVYDEEVEEPVLKAEAEKTEQEKTRQKEVDLIPKVQEAVNYGLQVLDSAFEQLDIKAG NSDSEEDDANGRVELILEPKDLYIDRPLPYLIGSKLFME
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Concentration	0.2 mg/ml
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C.

BACKGROUND

Keywords	FAM21A; family with sequence similarity 21, member A; FAM21B; bA98I6.1; bA56A21.1; WASH complex subunit FAM21A; WASH complex subunit FAM21B; family with sequence
-----------------	--

similarity 21, member B;

GENE INFORMATION

Entrez Gene ID [387680](#)

UniProt ID [Q641Q2](#)
