



Rabbit anti-Human FAM228B Polyclonal antibody (CABT-B9621)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: FDLKPLARAPYLLESQEEETVIYKNKGSSFLEREPLCYQEGNNPSAKEAISEGYFSSLS LSRNGRKTRMGLP
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Concentration	0.2 mg/ml
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C.

BACKGROUND

Keywords	FAM228B; family with sequence similarity 228, member B; protein FAM228B; UPF0638 protein B;
----------	---

GENE INFORMATION

Entrez Gene ID [375190](#)

UniProt ID [A0A087WZA1](#)
