



Rabbit anti-Human DESI1 Polyclonal antibody (CABT-B9623)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: IMLGKQLEGIWHTSIVVHKDEFFFGSGGISSCPPGGTLLGPPDSVVDVGSTEVTEE
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Concentration	0.1 mg/ml
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C.

BACKGROUND

Keywords	DESI1; desumoylating isopeptidase 1; DESI2; DeSI-1; PPPDE2; FAM152B; D15Wsu75e; DJ347H13.4; desumoylating isopeptidase 2; PPPDE peptidase domain containing 2; PPPDE peptidase domain-containing protein 2; family with sequence similarity 152, member B;
-----------------	--

GENE INFORMATION

Entrez Gene ID [27351](#)

UniProt ID [Q6ICB0](#)
