



Rabbit anti-Human FAM179B Polyclonal antibody (CABT-B9624)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: LLEHLKHKHSRVREEVVICISLLTYPSEDFDLPKLSFDLAPALVDSKRRVRQAALFAVLASSMGSGKTSILFKAVDTVELQDNGDGVMAVQ
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Concentration	0.05 mg/ml
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative.
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C.

BACKGROUND

Keywords	FAM179B; family with sequence similarity 179, member B; KIAA0423; protein FAM179B;
-----------------	--

GENE INFORMATION

Entrez Gene ID [23116](#)

UniProt ID [Q6P183](#)
