



Rabbit anti-Human CC SER1 Polyclonal antibody (CABT-B9626)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: S E R E E N I G L Q N G E T M L G T N S P R K L G F Y E Q H K A I A E H V K G I H P I S D S K I I P T S G D H H I F N K T S H G Y E A N P A K V L A S S L S P F R E G R F I
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	IHC, WB
Format	Liquid
Concentration	0.3 mg/ml
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2). 0.02% sodium azide is added as preservative
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C.

BACKGROUND

Keywords	CC SER1; coiled-coil serine-rich protein 1; FAM190A; serine-rich coiled-coil domain-containing protein 1; protein FAM190A; family with sequence similarity 190, member A;
-----------------	---

GENE INFORMATION

Entrez Gene ID [401145](#)

UniProt ID [Q9C0I3](#)
